EasyRoller

**100% Thermoplastique**

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	<b>TEST REPORT</b>	Doc. Nr.: 100127651 Date: July 26 <sup>th</sup> , 2010 Page 1 of 2	
<b>Document Number:</b> 100127651			
<b>Title:</b> Compatibility test between CEIA Walk-Through Metal Detectors and EasyRoller2			
<b>Description of the test</b>			
WTMD models used for the test	CEIA 02PN20 WTMD for airport application CEIA SMD601 WTMD for very high sensitivity application		
Item under test	EasyRoller2		
Description of the test	<p>Standard wheelchair manufactured with metal components are usually generating an alarm on the Walk-Through Metal Detector during the transit therefore the person on the wheel chair cannot be properly screened using the Walk-Through Metal Detector.</p> <p>EasyRoller 2 is a special wheelchair completely metal free according to the manufacturer specifications.</p> <p>The goal of this test has been to verify the compatibility of the CEIA Walk-Through Metal Detectors models 02PN20 and SMD601 with the EasyRoller2.</p> <p>In order to verify the signals measured by the Walk-Through Metal Detector, a laptop with CEIA MD-SCOPE application has been used.</p> <p>The laptop was connected to the Walk-Through Metal Detector via RS-232.</p>		
	 <small>CEIA 02PN20 for airport application</small>	 <small>CEIA SMD601 for high sensitivity application</small>	 <small>EasyRoller2</small>
<small>This document is property of CEIA S.p.A.. All right reserved</small>			
<small>MOD. AQ0020rev1</small>			

	<b>TEST REPORT</b>	Doc. Nr.: 100127651 Date: July 26 <sup>th</sup> , 2010 Page 2 of 2
	<p>The two CEIA Metal Detectors has been set at the following level of sensitivity:</p> <p>02PN20 – European Airport Security Level SMD601 – NU0601.02 Standard – Small Object Size Class (high sensitivity application)</p> <p>The EasyRoller2 was passed through both Metal Detectors several times verifying that for each transit no signal was measured by the Metal Detector.</p>	
Test results and final notes	<p>The result of the compatibility test between CEIA Walk-Through Metal Detectors and EasyRoller2 has been successful.</p> <p>It has been verified that no signal is generated in the Walk-Through Metal Detectors during the transit of the Easyroller2.</p> <p>This means that the Easyroller2 can be adopted to carry out security inspections of people who are not able to pass through the Walk-Through Metal Detector without the aid of a wheelchair.</p>	
Equipment, verification tools and instruments used	<ul style="list-style-type: none"> <li>• CEIA 02PN20 for airport application</li> <li>• CEIA SMD601 for high sensitivity application</li> <li>• EasyRoller2</li> <li>• CEIA MDSCOPE software</li> <li>• RS-232 serial cable</li> <li>• Laptop</li> </ul>	
	 <small>Sales Department Ing. Luca Pini</small>	
<small>This document is property of CEIA S.p.A.. All right reserved</small>		
<small>MOD. AQ0020rev1</small>		